

## Customer Information Bulletin

Part Number	Product Description	Issued by	Hellie King, Product Manager
120-112XSFV	1.2micron filter microset with manual priming valve	Document No	CB17-008
		Date	24 April 2017

### Purpose of Bulletin

To communicate a product enhancement to the above giving set and a product code change from **120-112XSFV** to **120-112XSFVK**.

### Detail

In order to improve the ease and ensure the correct placement of the set, into the BodyGuard pump, there has been an addition of a circular blue key.


PN120-112XSFV currently has the below segment with one key:



The image below shows the addition, on PN120-112XSFVK, of a blue circular key, referred to as the '2<sup>nd</sup> key'.



The image below is a representation of how the segment is loaded into the pump. The blue circular key is loaded first followed by the black key.



If you have any further questions or enquiries, please contact Customer Experience Team on 01253 206700 or email [customersupport@cmemmedical.co.uk](mailto:customersupport@cmemmedical.co.uk)